

HYPERTENSION

Use the TONIC SUPER-IMPOSED WAVE of the No. 1510 Modernistic Sinustat.

TECHNIQUE: Patient sitting on treatment table ready to lie down. Place one 3" x 4" Sinusoidal electrode, thoroughly moistened and soaped, directly over the third and fourth dorsal vertebrae. The other 3" x 4" electrode directly over the 10th, 11th and 12th dorsal vertebrae. Back up the electrodes with rubber sponges and have patient lie on back, making firm contact with both electrodes. Gradually turn on the No. 4 Tonic Super-Imposed Wave of the Sinustat to patient's comfortable tolerance. Use a pulsation rate of about 20 to 40 per minute. Treat about 15 to 20 minutes. Repeat every other day or daily if condition is serious or slow to respond. Ten to twelve treatments usually produce results but continue until relief seems permanent.

It is a good plan to take the blood pressure before and after treatment and note change.

CAUTION: Do not use this treatment for hypertension, if the pulse pressure is below 20 millimeters.

By this technique, we stimulate the function of the centers of inhibition and dilatation, producing a general systemic effect.